

E 900-6
Mobile unit
for the optimum cleaning of instruments on 6 levels.



- Equipped with connections for hollow instruments
- Loading capacity: 18 DIN mesh trays*
- Equipped with 3 integrated spray arms
- The loading areas for the individual levels are different
- Magnetic strip for mobile unit recognition (incl. magnets ML)

* DIN mesh tray 485 x 254 x 50 mm

EAN: 4002515263358 / Material number: 09525120 / Old Material Number: 69590007D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
Mobile unit for MIS instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
Capacity	
DIN mesh trays [number]	18
SPRI mesh trays [number]	12
Lumened instruments [number]	16
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	553
External dimensions, net width in mm	636
External dimensions, net depth in mm	790
External dimensions, gross height in mm	680
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	650
Net weight in kg	12.85
Gross weight in kg	12.85
Companion pack	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1