

A 803  
Holder  
for the optimum loading of transmission instruments.

EAN: 4002515460993 / Material number: 09863100 / Old Material Number: 69680300D



- For use in A 105/1 upper basket along with A 800 tubular filter
- Suitable for numerous contra-angles & hand pieces as well as turbines

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
<b>Product category</b>	
Holder	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of dental instruments	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	32
External dimensions, net width in mm	32
External dimensions, net depth in mm	40
External dimensions, gross height in mm	80
External dimensions, gross width in mm	95
External dimensions, gross depth in mm	110
Net weight in kg	0.05
Gross weight in kg	0.09