

UG 30-60/60-85

Box plinth

to raise the machine, 60 cm deep, 60 cm wide.

EAN: 4002515558331 / Material number: 10087350 / Old Material Number: 69306001D



- Suitable for all machines with a width and depth of 60 cm

Product category	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8504	•
PG 8583	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Plinth/base	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in mm	317
External dimensions, net width in mm	600
External dimensions, net depth in mm	600
External dimensions, gross height in mm	310
External dimensions, gross width in mm	610
External dimensions, gross depth in mm	610
Net weight in kg	9.70
Gross weight in kg	11.10