

UBS 3

Conversion kit

For use with selected series G 78 upper baskets in PG 85 machines.

EAN: 4002515588994 / Material number: 10131110 / Old Material Number: 69215016D



- For older load carriers without dryer connection
- For older load carriers with welded water supply

Product category	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Conversion kit	•
Product properties	
Colour	Stainless Steel, Black
Dimensions and weight	
External dimensions, net height in mm	90
External dimensions, net width in mm	305
External dimensions, net depth in mm	176
External dimensions, gross height in mm	100
External dimensions, gross width in mm	310
External dimensions, gross depth in mm	210
Net weight in kg	0.59
Gross weight in kg	0.69