

DE-CS7-85

Appliance lid

for the top cover on the machine, 70 cm deep, 90 cm wide.



EAN: 4002515648155 / Material number: 10222780 / Old Material Number: 69216011D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8582 CD	•
PG 8583 CD	•
Product category	
Lid	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in mm	15
External dimensions, net width in mm	700
External dimensions, net depth in mm	900
External dimensions, gross height in mm	80
External dimensions, gross width in mm	790
External dimensions, gross depth in mm	990
Net weight in kg	14.45
Gross weight in kg	16.75