

Set GF-1 100 mm  
 Set of machine feet  
 Can be adjusted to the optimum operating height.



- Suitable for all dishwashers in the PG 8172 range
- Comprising 4 machine feet
- To raise by 10 cm

EAN: 4002515641866 / Material number: 10249450 / Old Material Number: 67820030D

<b>Product category</b>	
Throughfeed tank dishwashers	•
<b>Machine affiliation</b>	
PG 8172	•
PG 8172 Eco	•
PTD 901	•
<b>Product category</b>	
Other accessories	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	90
External dimensions, gross width in mm	150
External dimensions, gross depth in mm	280
Net weight in kg	2.40
Gross weight in kg	2.70