

U 403
Basket
for the optimum cleaning of cutlery.



- Close-meshed base
- Also suitable for sundae dishes
- Without cutlery inserts

EAN: 4002515691175 / Material number: 10358590 / Old Material Number: 67640301D

Product category	
Tank dishwashers	•
Machine affiliation	
PG 8164	•
PTD 701	•
Product category	
Basket for cutlery	•
Basket for saucepans/bowls, etc.	•
Product properties	
Material	Wire, rilsan-coated
Colour	White
Dimensions and weight	
External dimensions, net height in mm	125
External dimensions, net width in mm	400
External dimensions, net depth in mm	400
External dimensions, gross height in mm	170
External dimensions, gross width in mm	520
External dimensions, gross depth in mm	620
Net weight in kg	2.00
Gross weight in kg	3.00