

U 404

Basket

For the optimum loading of cutlery.



- With 6 plastic cutlery holders
- Compartment size 106 x 106 mm

EAN: 4002515691182 / Material number: 10358610 / Old Material Number: 67640401D

<b>Product category</b>	
Tank dishwashers	•
<b>Machine affiliation</b>	
PG 8164	•
PTD 701	•
<b>Product category</b>	
Basket for cutlery	•
<b>Product properties</b>	
Material	Wire, rilsan-coated, Plastic
Colour	Grey, White
Number of compartments	6
Compartment width in mm	106
Compartment length in mm	106
<b>Dimensions and weight</b>	
External dimensions, net height in mm	112
External dimensions, net width in mm	400
External dimensions, net depth in mm	400
External dimensions, gross height in mm	170
External dimensions, gross width in mm	520
External dimensions, gross depth in mm	620
Net weight in kg	1.65
Gross weight in kg	2.80