

A 815 Nozzle

Optimum reprocessing of transmission instruments with external spray channel.

EAN: 4002515830413 / Material number: 10619690 / Old Material Number: 68681500D



- For use with A 812, A 813 or in A 105/1 upper basket
- Consists of 4 nozzles with silicone hoses, L approx. 114 mm

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Transparent
Area of application	
Reprocessing of dental instruments	•
Capacity	
Transmission instruments [number]	4
Dimensions and weight	
External dimensions, net height in mm	10
External dimensions, net width in mm	10
External dimensions, net depth in mm	114
External dimensions, gross height in mm	80
External dimensions, gross width in mm	95
External dimensions, gross depth in mm	110
Diameter in mm	10
Net weight in kg	0.05
Gross weight in kg	0.12