

A 833

Injector nozzle

for the optimum cleaning of a variety of load items.



- Suitable for mobile injector units, with screw connection
- For optimum cleaning results thanks to lateral drill holes

EAN: 4002515941225 / Material number: 10635480 / Old Material Number: 69745302D

Product classification	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	4
External dimensions, net depth in mm	110
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.00
Gross weight in kg	0.01