

A 834

Injector nozzle

for the optimum cleaning of diverse loads, particularly for proctoscopes.

EAN: 4002515940709 / Material number: 10635650 / Old Material Number: 69745403D



- Suitable for proctoscopes, trocar sleeves and other load items
- Equipped with a spring for height adjustment
- For optimum cleaning results thanks to lateral drill holes

Product classification	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	8
External dimensions, net depth in mm	152
External dimensions, gross height in mm	12
External dimensions, gross width in mm	120
External dimensions, gross depth in mm	170
Net weight in kg	0.03
Gross weight in kg	0.03