

A 606
Module
for the use of e.g. lab glassware inserts on the top level of A 503 mobile units



- With spray arm

EAN: 4002516050537 / Material number: 10973990 / Old Material Number: 69660601D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Other inserts/other modules	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	97
External dimensions, net width in mm	505
External dimensions, net depth in mm	606
External dimensions, gross height in mm	120
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	620
Net weight in kg	2.40
Gross weight in kg	3.00