


A 605
Module
for optimum reprocessing of vials and test tubes, etc.



- Equipped with 98 injector nozzles
- Max. H 160 mm, max. Ø 20 mm
- For use in A 503 mobile unit

EAN: 4002516063339 / Material number: 10999010 / Old Material Number: 69660500D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert/module for test tubes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of pipettes/phials	•
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	98
Phials [number]	98
Dimensions and weight	
External dimensions, net height in mm	158
External dimensions, net width in mm	254
External dimensions, net depth in mm	604
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	4.40
Gross weight in kg	6.40