

A 612  
Module  
for the optimum reprocessing of large-volume laboratory glassware.



- Can be equipped with injector nozzles based on individual needs
- Suitable for 2–10 l laboratory flasks, 2–6 l round flasks
- Suitable for 2–5 l Erlenmeyer flasks, 2–5 l measuring flasks
- For use in A 503 mobile unit

EAN: 4002516063353 / Material number: 10999050 / Old Material Number: 69661200D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Insert/module for narrow-necked glassware	•
Module for large-volume laboratory glassware	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Large-volume laboratory glassware [number]	4
<b>Dimensions and weight</b>	
External dimensions, net height in mm	221
External dimensions, net width in mm	253
External dimensions, net depth in mm	608
External dimensions, gross height in mm	245
External dimensions, gross width in mm	270
External dimensions, gross depth in mm	635
Net weight in kg	3.60
Gross weight in kg	4.40