

A 506  
Mobile unit  
with 4 levels for cleaning laboratory glassware and utensils.



- For holding inserts and mesh trays

EAN: 4002516063223 / Material number: 10999130 / Old Material Number: 69650600D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Mobile unit for inserts	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
<b>Capacity</b>	
DIN mesh trays [number]	12
SPRI mesh trays [number]	8
<b>Dimensions and weight</b>	
External dimensions, net height in mm	555
External dimensions, net width in mm	650
External dimensions, net depth in mm	790
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	4.40
Gross weight in kg	6.40