


APWD 060

Basket

for the optimum loading of different inserts and mesh trays.



- Vertical clearance 230 mm, width 370 mm, depth 395 mm

EAN: 4002516066644 / Material number: 11017140 / Old Material Number: 68406000D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PWD 8531	•
PWD 8532	•
PWD 8534	•
Product category	
Other basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, White
Area of application	
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of orthopaedic surgery instruments	•
Reprocessing of trays	•
Dimensions and weight	
External dimensions, net height in mm	85
External dimensions, net width in mm	395
External dimensions, net depth in mm	469
External dimensions, gross height in mm	262
External dimensions, gross width in mm	408
External dimensions, gross depth in mm	512
Net weight in kg	1.42
Gross weight in kg	2.92