


APWD 061

Basket

for the optimum loading of handpieces, different inserts and mesh trays.

EAN: 4002516066651 / Material number: 11017150 / Old Material Number: 68406100D



- Equipped with 6 holders for dental/ surgical handpieces & turbines
- Equipped with 6 injectors Ø 3.8 x 15 mm for extractors
- Active inner drying of dental and surgical handpieces
- Loading dimensions H 230, W 320, D 395 mm

Product category	
Thermal disinfectors	•
Machine affiliation	
PWD 8531	•
PWD 8532	•
Product category	
Injector basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
Area of application	
Reprocessing of dental instruments	•
Capacity	
Transmission instruments [number]	6
Lumened instruments [number]	6
Dimensions and weight	
External dimensions, net height in mm	85
External dimensions, net width in mm	395
External dimensions, net depth in mm	469
External dimensions, gross height in mm	262
External dimensions, gross width in mm	408
External dimensions, gross depth in mm	512
Net weight in kg	2.18
Gross weight in kg	3.68