

APWD 062

Insert

for the optimum loading of upright instruments.



- Suitable for probes, tweezers and mouth specula
- Suitable for loading files, chisels and podiatry instruments.

EAN: 4002516066668 / Material number: 11017160 / Old Material Number: 68406200D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PWD 8531	•
PWD 8532	•
PWD 8534	•
PWD 8682 C	•
Product category	
Insert for instruments	•
Insert/module for dental instruments	•
Area of application	
Reprocessing of dental instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in mm	122
External dimensions, net width in mm	124
External dimensions, net depth in mm	395
External dimensions, gross height in mm	150
External dimensions, gross width in mm	140
External dimensions, gross depth in mm	425
Net weight in kg	0.71
Gross weight in kg	1.00