

APWD 063

Insert

for the optimum loading of up to 8 trays or cassettes.

EAN: 4002516066675 / Material number: 11017170 / Old Material Number: 68406300D



- Equipped with 9 holders, spacing: 40 mm

| | |
|---|-----------------|
| Product category | |
| Thermal disinfectors | • |
| Machine affiliation | |
| PWD 8531 | • |
| PWD 8532 | • |
| Product category | |
| Insert for trays | • |
| Insert/module for dental instruments | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of trays | • |
| Capacity | |
| Trays [quantity] | 8 |
| Dimensions and weight | |
| External dimensions, net height in mm | 110 |
| External dimensions, net width in mm | 103 |
| External dimensions, net depth in mm | 389 |
| External dimensions, gross height in mm | 185 |
| External dimensions, gross width in mm | 222 |
| External dimensions, gross depth in mm | 556 |
| Net weight in kg | 0.30 |
| Gross weight in kg | 1.09 |