


APWD 067

for optimum reprocessing of transmission instruments in APWD 061 basket

for optimum reprocessing of transmission instruments in APWD 061 basket



- Can only be used with matching APWD 068 lower part

EAN: 4002516066729 / Material number: 11017200 / Old Material Number: 68406700D

Product category	
Thermal disinfectors	•
Machine affiliation	
PWD 8531	•
PWD 8532	•
Product category	
Adapter	•
Holder and/or adapter for transmission instruments	•
Product properties	
Material	Plastic
Colour	Blue
Area of application	
Reprocessing of dental instruments	•
Dimensions and weight	
External dimensions, net height in mm	2
External dimensions, net width in mm	346
External dimensions, net depth in mm	3
External dimensions, gross height in mm	12
External dimensions, gross width in mm	120
External dimensions, gross depth in mm	170
Net weight in kg	0.01
Gross weight in kg	0.01