

A 623
Module
for the optimum loading of measuring cylinders.



- Equipped with 4 nozzles
- For measuring cylinders with 1–2 l capacity
- Pivoting supports for ease of use

EAN: 4002516071440 / Material number: 11025560 / Old Material Number: 69662300D

| | |
|--|--------------------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |
| Product category | |
| Insert/module for measuring cylinder | • |
| Product properties | |
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel, Black |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Reprocessing of measuring cylinders | • |
| Capacity | |
| Measuring cylinders [number] | 4 |
| Dimensions and weight | |
| External dimensions, net height in mm | 437 |
| External dimensions, net width in mm | 255 |
| External dimensions, net depth in mm | 608 |
| External dimensions, gross height in mm | 540 |
| External dimensions, gross width in mm | 530 |
| External dimensions, gross depth in mm | 500 |
| Net weight in kg | 6.30 |
| Gross weight in kg | 7.80 |