

A 838

Key set

For adjusting basket heights and installing/removing injector nozzles.



- Set consisting of Torx wrench and open spanner
- Suitable for Miele upper baskets and mobile units

EAN: 4002516171720 / Material number: 11054290 / Old Material Number: 69783800D

Product classification	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8505	•
Product category	
Other accessories	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.00
Gross weight in kg	0.00