

A 843

Injector nozzle

with plastic supports, nozzle diameter 4 mm, length 185 mm, 1 item



- Suitable for mobile injector units & baskets screw connection
- Part of the EasyLoad system
- Packaging unit with 1 nozzle

EAN: 4002516188773 / Material number: 11056440 / Old Material Number: 69684300D

| | |
|--|--------------------------|
| Product category | |
| Washer-disinfectors, laboratory | • |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| G 7825 | • |
| PG 8527 | • |
| PG 8528 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8504 | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8593 | • |
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |
| Product category | |
| Injector nozzle | • |
| Product properties | |
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel, Black |
| Part of the EasyLoad system | • |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Dimensions and weight | |
| External dimensions, net height in mm | 185 |
| External dimensions, net width in mm | 76 |
| External dimensions, net depth in mm | 71 |
| External dimensions, gross height in mm | 12 |
| External dimensions, gross width in mm | 80 |
| External dimensions, gross depth in mm | 120 |
| Net weight in kg | 0.02 |
| Gross weight in kg | 0.02 |