


A 862
Holder rack
for injector modules A 302/2 and A 302/3



- Arrangement: 4 x 8 fields
- For use on A 302/2 and A 302/3 injector modules
- For the optimum positioning of laboratory glassware
- Plastic-coated

EAN: 4002516087687 / Material number: 11056820 / Old Material Number: 69686200D

| | |
|--|---|
| Product category | |
| Washer-disinfectors, laboratory | • |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PG 8504 | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8593 | • |
| Product category | |
| Other inserts/other modules | • |
| Product properties | |
| Material | Stainless steel, Plastic, Wire, rilsan-coated |
| Colour | Stainless Steel, Black |
| Number of compartments | 32 |
| Compartment width in mm | 56 |
| Compartment length in mm | 52 |
| Part of the EasyLoad system | • |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0-0 |
| Frequency in Hz | 0-0 |
| Total rated load in kW | 0.00-0.00 |
| Fuse rating in A | 0-0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 30 |
| External dimensions, net width in mm | 250 |
| External dimensions, net depth in mm | 350 |
| External dimensions, gross height in mm | 35 |
| External dimensions, gross width in mm | 300 |
| External dimensions, gross depth in mm | 400 |
| Net weight in kg | 0.10 |
| Gross weight in kg | 0.18 |