


APLW 000

Mobile unit

for optimum loading of laboratory glassware on upper levels.

EAN: 4002516115281 / Material number: 11104470 / Old Material Number: 69300000D



- Equipped with 1 integrated spray arm
- For holding inserts
- Clearance: 484 x 525 mm

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Insert/module for beakers	•
Upper basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	160
External dimensions, net width in mm	560
External dimensions, net depth in mm	532
External dimensions, gross height in mm	190
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	610
Net weight in kg	4.30
Gross weight in kg	6.07