


APLW 002

Mobile unit

for optimum loading of large-volume laboratory glassware, 5 and 10 l.

EAN: 4002516115304 / Material number: 11104490 / Old Material Number: 69300200D



- Equipped with 5 nozzles
- For holding up to 4 flasks, max. Ø 240 mm
- For holding up to 5 flasks, max. Ø 190 mm

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Module for large-volume laboratory glassware	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Flasks [number]	5
<b>Dimensions and weight</b>	
External dimensions, net height in mm	413
External dimensions, net width in mm	535
External dimensions, net depth in mm	582
External dimensions, gross height in mm	580
External dimensions, gross width in mm	565
External dimensions, gross depth in mm	615
Net weight in kg	5.94
Gross weight in kg	8.49