

APLW 003

Mobile unit

for optimum loading of large-volume laboratory glassware, up to 20 l.



- Equipped with 4 nozzles
- For holding up to 2 flasks, max. Ø 280 mm
- For loading up to 2 bottles with a max. Ø of 190 mm

EAN: 4002516115410 / Material number: 11104500 / Old Material Number: 69300300D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Flasks [number]	2
Erlenmeyer flasks [number]	2
Dimensions and weight	
External dimensions, net height in mm	400
External dimensions, net width in mm	535
External dimensions, net depth in mm	582
External dimensions, gross height in mm	580
External dimensions, gross width in mm	565
External dimensions, gross depth in mm	615
Net weight in kg	5.08
Gross weight in kg	7.63