

APLW 004

Mobile unit

for optimum loading of large-volume laboratory glassware, up to 50 l.

EAN: 4002516115427 / Material number: 11104510 / Old Material Number: 69300400D



- Equipped with 1 nozzle

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Module for large-volume laboratory glassware	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Flasks [number]	1
<b>Dimensions and weight</b>	
External dimensions, net height in mm	588
External dimensions, net width in mm	535
External dimensions, net depth in mm	582
External dimensions, gross height in mm	760
External dimensions, gross width in mm	625
External dimensions, gross depth in mm	715
Net weight in kg	4.68
Gross weight in kg	9.83