

APLW 005
 Mobile unit
 for optimum loading of up to 56 pipettes.



- For holding up to 56 pipettes arranged diagonally

EAN: 4002516115434 / Material number: 11104520 / Old Material Number: 69300500D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Mobile unit for pipettes	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel
Area of application	
Reprocessing of pipettes/phials	•
Capacity	
Pipettes [number]	56
Dimensions and weight	
External dimensions, net height in mm	308
External dimensions, net width in mm	535
External dimensions, net depth in mm	582
External dimensions, gross height in mm	570
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	610
Net weight in kg	6.40
Gross weight in kg	9.20