

APLW 010
Upper basket
with 156 unoccupied injector slots for nozzles with Ø 2.5 mm.



- Only suitable for use on upper levels

EAN: 4002516115489 / Material number: 11104570 / Old Material Number: 69301000D

| | |
|--|------------------------------|
| Product classification | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Injector basket | • |
| Upper basket | • |
| Product properties | |
| Material | Stainless steel, Plastic |
| Colour | Blue, Stainless Steel, White |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Narrow-necked glassware [number] | 156 |
| Dimensions and weight | |
| External dimensions, net height in mm | 162 |
| External dimensions, net width in mm | 534 |
| External dimensions, net depth in mm | 569 |
| External dimensions, gross height in mm | 355 |
| External dimensions, gross width in mm | 565 |
| External dimensions, gross depth in mm | 615 |
| Net weight in kg | 5.16 |
| Gross weight in kg | 7.58 |