

APLW 013/1  
Upper basket  
equipped with 42 APLW 051 injector nozzles



**Image  
in progress**

- Only suitable for use on upper levels

EAN: 4002516115526 / Material number: 11104610 / Old Material Number: 69301301D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-necked glassware [number]	42
Flasks [number]	42
Erlenmeyer flasks [number]	42
Round flasks [number]	42
Dimensions and weight	
External dimensions, net height in mm	160
External dimensions, net width in mm	534
External dimensions, net depth in mm	569
External dimensions, gross height in mm	355
External dimensions, gross width in mm	565
External dimensions, gross depth in mm	615
Net weight in kg	5.88
Gross weight in kg	8.28
Companion pack	
APLW 051, injector nozzle Ø 4 mm, length 110 mm [number]	42