

APLW 015  
Upper basket  
with 16 unoccupied injector slots.



**Image  
in progress**

- Only suitable for use on upper levels

EAN: 4002516115571 / Material number: 11104650 / Old Material Number: 69301500D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-necked glassware [number]	16
Flasks [number]	16
Erlenmeyer flasks [number]	16
Round flasks [number]	16
<b>Dimensions and weight</b>	
External dimensions, net height in mm	160
External dimensions, net width in mm	532
External dimensions, net depth in mm	560
External dimensions, gross height in mm	340
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	600
Net weight in kg	4.00
Gross weight in kg	6.60