

APLW 018
Upper basket
with 121 + 24 unoccupied injector slots.



**Image
in progress**

- Only suitable for use on upper levels

EAN: 4002516115601 / Material number: 11104680 / Old Material Number: 69301800D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	121
Phials [number]	121
Narrow-necked glassware [number]	24
Flasks [number]	24
Erlenmeyer flasks [number]	24
Round flasks [number]	24
Dimensions and weight	
External dimensions, net height in mm	162
External dimensions, net width in mm	534
External dimensions, net depth in mm	569
External dimensions, gross height in mm	355
External dimensions, gross width in mm	565
External dimensions, gross depth in mm	615
Net weight in kg	5.33
Gross weight in kg	7.76