


APLW 018/1

Upper basket

equipped with 121 APLW 048 injector nozzles and 24 APLW 051 injector nozzles.

EAN: 4002516115717 / Material number: 11104690 / Old Material Number: 69301801D



- Only suitable for use on upper levels

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	121
Phials [number]	121
Narrow-necked glassware [number]	24
Flasks [number]	24
Erlenmeyer flasks [number]	24
Round flasks [number]	24
<b>Dimensions and weight</b>	
External dimensions, net height in mm	200
External dimensions, net width in mm	535
External dimensions, net depth in mm	570
External dimensions, gross height in mm	355
External dimensions, gross width in mm	565
External dimensions, gross depth in mm	615
Net weight in kg	6.92
Gross weight in kg	11.42
<b>Companion pack</b>	
APLW 051, injector nozzle Ø 4 mm, length 110 mm [number]	32
APLW 048 injector nozzle, Ø 2.5 mm, length 80 mm [number]	121