

APLW 021

Lower basket

with 110 unoccupied injector slots for nozzles with Ø 2.5 mm.



**Image  
in progress**

- Only suitable for use on the lowest level

EAN: 4002516115762 / Material number: 11104740 / Old Material Number: 69302100D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	110
Phials [number]	110
Narrow-necked glassware [number]	110
Flasks [number]	110
Dimensions and weight	
External dimensions, net height in mm	575
External dimensions, net width in mm	504
External dimensions, net depth in mm	450
External dimensions, gross height in mm	355
External dimensions, gross width in mm	562
External dimensions, gross depth in mm	615
Net weight in kg	5.91
Gross weight in kg	8.35