

APLW 022

Lower basket

with 70 unoccupied injector slots for nozzles with Ø 2.5 mm.



**Image  
in progress**

- Only suitable for use on the lowest level

EAN: 4002516115779 / Material number: 11104750 / Old Material Number: 69302200D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-necked glassware [number]	70
Flasks [number]	70
Erlenmeyer flasks [number]	70
Round flasks [number]	70
<b>Dimensions and weight</b>	
External dimensions, net height in mm	580
External dimensions, net width in mm	515
External dimensions, net depth in mm	534
External dimensions, gross height in mm	355
External dimensions, gross width in mm	562
External dimensions, gross depth in mm	615
Net weight in kg	5.83
Gross weight in kg	8.27