

APLW 025  
Lower basket  
with 24 unoccupied injector slots for nozzles.



- Only suitable for use on the lowest level
- 230 x 490 mm usable space

EAN: 4002516115830 / Material number: 11104810 / Old Material Number: 69302500D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-necked glassware [number]	24
Flasks [number]	24
Erlenmeyer flasks [number]	24
Round flasks [number]	24
<b>Dimensions and weight</b>	
External dimensions, net height in mm	136
External dimensions, net width in mm	535
External dimensions, net depth in mm	582
External dimensions, gross height in mm	360
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	600
Net weight in kg	4.40
Gross weight in kg	7.40