

APLW 025/1

Lower basket

equipped with 12 APLW 056 injector nozzles and 12 APLW 052 injector nozzles.



**Image  
in progress**

- Only suitable for use on the lowest level
- 230 x 490 mm usable space

EAN: 4002516115847 / Material number: 11104820 / Old Material Number: 69302501D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-necked glassware [number]	24
Flasks [number]	24
Erlenmeyer flasks [number]	24
Round flasks [number]	24
<b>Dimensions and weight</b>	
External dimensions, net height in mm	270
External dimensions, net width in mm	535
External dimensions, net depth in mm	570
External dimensions, gross height in mm	350
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	610
Net weight in kg	5.25
Gross weight in kg	8.39
<b>Companion pack</b>	
APLW 052 injector nozzle, Ø 4 mm, length 175 mm [number]	12
APLW 056 injector nozzle, Ø 6 mm, length 225 mm [number]	12