

APLW 026  
Lower basket  
with 20 unoccupied injector slots for nozzles.



**Image  
in progress**

- Only suitable for use on the lowest level

EAN: 4002516115854 / Material number: 11104830 / Old Material Number: 69302600D

|  |                               |
|--|-------------------------------|
| <b>Product classification</b>                  |                               |
| Large-capacity washer-disinfectors, laboratory | •                             |
| <b>Machine affiliation</b>                     |                               |
| PLW 6011                                       | •                             |
| PLW 6111                                       | •                             |
| <b>Product category</b>                        |                               |
| Injector basket                                | •                             |
| Lower basket                                   | •                             |
| <b>Product properties</b>                      |                               |
| Material                                       | Stainless steel, Plastic      |
| Colour   | Stainless Steel, Black, White |
| <b>Area of application</b>                     |                               |
| Reprocessing of laboratory glassware           | •                             |
| <b>Capacity</b>                                |                               |
| Narrow-necked glassware [number]               | 20                            |
| Flasks [number]                                | 20                            |
| Erlenmeyer flasks [number]                     | 20                            |
| Round flasks [number]                          | 20                            |
| <b>Dimensions and weight</b>                   |                               |
| External dimensions, net height in mm          | 570                           |
| External dimensions, net width in mm           | 536                           |
| External dimensions, net depth in mm           | 516                           |
| External dimensions, gross height in mm        | 340                           |
| External dimensions, gross width in mm         | 555                           |
| External dimensions, gross depth in mm         | 615                           |
| Net weight in kg                               | 4.92                          |
| Gross weight in kg                             | 7.40                          |