


APLW 026/1
Lower basket
equipped with 20 APLW 056 injector nozzles.

EAN: 4002516115861 / Material number: 11104840 / Old Material Number: 69302601D



- Only suitable for use on the lowest level

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-necked glassware [number]	20
Flasks [number]	20
Erlenmeyer flasks [number]	20
Round flasks [number]	20
Dimensions and weight	
External dimensions, net height in mm	270
External dimensions, net width in mm	535
External dimensions, net depth in mm	570
External dimensions, gross height in mm	340
External dimensions, gross width in mm	555
External dimensions, gross depth in mm	615
Net weight in kg	4.92
Gross weight in kg	9.20
Companion pack	
APLW 056 injector nozzle, Ø 6 mm, length 225 mm [number]	20