

APLW 032
 Lower basket
 with 121 + 24 unoccupied injector slots for nozzles.



- Only suitable for use on the lowest level

EAN: 4002516115939 / Material number: 11104910 / Old Material Number: 69303200D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	121
Phials [number]	121
Narrow-necked glassware [number]	24
Flasks [number]	24
Dimensions and weight	
External dimensions, net height in mm	570
External dimensions, net width in mm	536
External dimensions, net depth in mm	70
External dimensions, gross height in mm	355
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	610
Net weight in kg	5.99
Gross weight in kg	8.31