

APLW 032/1

Lower basket

equipped with 121 APLW 048 injector nozzles and 24 APLW 056 injector nozzles.



**Image  
in progress**

- Only suitable for use on the lowest level

EAN: 4002516115946 / Material number: 11104920 / Old Material Number: 69303201D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	121
Phials [number]	121
Narrow-necked glassware [number]	24
Flasks [number]	24
<b>Dimensions and weight</b>	
External dimensions, net height in mm	270
External dimensions, net width in mm	535
External dimensions, net depth in mm	575
External dimensions, gross height in mm	355
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	610
Net weight in kg	5.99
Gross weight in kg	13.98
<b>Companion pack</b>	
APLW 048 injector nozzle, Ø 2.5 mm, length 80 mm [number]	121
APLW 056 injector nozzle, Ø 6 mm, length 225 mm [number]	24