

APLW 038

Insert

for the optimum loading of wide-necked Erlenmeyer flasks & measuring cylinders.



- Equipped with 14 spring hooks, 175 mm & 14 105 mm height

EAN: 4002516116004 / Material number: 11104980 / Old Material Number: 69303800D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Insert/module for narrow-necked glassware	•
Insert for various flasks	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	185
External dimensions, net width in mm	195
External dimensions, net depth in mm	450
External dimensions, gross height in mm	345
External dimensions, gross width in mm	435
External dimensions, gross depth in mm	510
Net weight in kg	0.60
Gross weight in kg	2.41