

APLW 039

Insert

for optimum loading of half petri dishes.



- Suitable for 26 half petri dishes
- Spacing between the holders 24 mm

EAN: 4002516116011 / Material number: 11104990 / Old Material Number: 69303900D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Insert/module for half petri dishes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	83
External dimensions, net width in mm	237
External dimensions, net depth in mm	430
External dimensions, gross height in mm	255
External dimensions, gross width in mm	410
External dimensions, gross depth in mm	510
Net weight in kg	1.17
Gross weight in kg	2.47