


APLW 040

Insert

for optimum loading of test tubes, height 100 mm.

EAN: 4002516116028 / Material number: 11105000 / Old Material Number: 69304000D



- For use with APLW 043 lid

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Insert/module for test tubes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Maximum vertical clearance in mm	115
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	100
External dimensions, net width in mm	220
External dimensions, net depth in mm	220
External dimensions, gross height in mm	330
External dimensions, gross width in mm	315
External dimensions, gross depth in mm	315
Net weight in kg	0.47
Gross weight in kg	1.58