


APLW 042

Insert

for optimum loading of test tubes, height 200 mm.



- For use with APLW 043 lid

EAN: 4002516116042 / Material number: 11105020 / Old Material Number: 69304200D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Insert/module for test tubes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Maximum vertical clearance in mm	215
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	202
External dimensions, net width in mm	220
External dimensions, net depth in mm	220
External dimensions, gross height in mm	325
External dimensions, gross width in mm	305
External dimensions, gross depth in mm	310
Net weight in kg	0.65
Gross weight in kg	1.60