

APLW 042

Insert

for optimum loading of test tubes, height 200 mm.

EAN: 4002516116042 / Material number: 11105020 / Old Material Number: 69304200D



- For use with APLW 043 lid

|  |                 |
|--|-----------------|
| <b>Product classification</b>                  |                 |
| Large-capacity washer-disinfectors, laboratory | •               |
| <b>Machine affiliation</b>                     |                 |
| PLW 6011                                       | •               |
| PLW 6111                                       | •               |
| <b>Product category</b>                        |                 |
| Insert/module for test tubes                   | •               |
| <b>Product properties</b>                      |                 |
| Material                                       | Stainless steel |
| Colour   | Stainless Steel |
| Maximum vertical clearance in mm               | 215             |
| <b>Area of application</b>                     |                 |
| Reprocessing of laboratory glassware           | •               |
| <b>Dimensions and weight</b>                   |                 |
| External dimensions, net height in mm          | 202             |
| External dimensions, net width in mm           | 220             |
| External dimensions, net depth in mm           | 220             |
| External dimensions, gross height in mm        | 325             |
| External dimensions, gross width in mm         | 305             |
| External dimensions, gross depth in mm         | 310             |
| Net weight in kg                               | 0.65            |
| Gross weight in kg                             | 1.60            |