


APLW 044
 Separation net
 for dividing the base area in APLW 040–042.



EAN: 4002516116066 / Material number: 11105040 / Old Material Number: 69304400D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Insert/module for test tubes	•
Holder	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	75
External dimensions, net width in mm	205
External dimensions, net depth in mm	200
External dimensions, gross height in mm	215
External dimensions, gross width in mm	170
External dimensions, gross depth in mm	270
Net weight in kg	0.36
Gross weight in kg	0.60