


APLW 044  
 Separation net  
 for dividing the base area in APLW 040–042.



EAN: 4002516116066 / Material number: 11105040 / Old Material Number: 69304400D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Insert/module for test tubes	•
Holder	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	75
External dimensions, net width in mm	205
External dimensions, net depth in mm	200
External dimensions, gross height in mm	215
External dimensions, gross width in mm	170
External dimensions, gross depth in mm	270
Net weight in kg	0.36
Gross weight in kg	0.60