

APLW 048

Injector nozzle

with plastic support, star-shaped, Ø 2.5, length 80 mm.

EAN: 4002516116097 / Material number: 11105080 / Old Material Number: 69304800D



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 4 mm

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	32
External dimensions, net width in mm	32
External dimensions, net depth in mm	90
External dimensions, gross height in mm	5
External dimensions, gross width in mm	110
External dimensions, gross depth in mm	167
Net weight in kg	0.00
Gross weight in kg	0.00